

- **Consumables**

- *Infection control kit, XO FLOW, FL-464, see details below*
- *Dashboard protection foil, MN-880, see details below*
- *Air instrument return air oil separator sponge, 10 pcs, AP-983*
- *XO FLOW Service kit, year 2, 4 etc., AR-132*
- *XO Water softener, filter cartridge, UH-200*
- *High volume and saliva suction hoses complete (175 cm) for right-handed operator, AR-124*
- *High volume and saliva suction hoses complete (260 cm) for right- and left-handed operator, AR-127*
- *Suction filter with filter holder, 20 pcs, AR-209*
- *Handles, XO FLEX + XO FLOW, 20 pcs, AR-299*
- *Peroxide test strips, 100 pcs, UH-238*
- *Silicone grease, O-rings (NBR & EPDM), 25 ml, UG-928*
- *Darina R2 grease, 100 ml, YR-026*
- *Paraliq grease, 100 ml, YR-032*
- *Sunaplex 781 grease, 1000 ml, FA-385*
- *Sunaplex 781 grease, 100 ml, YR-033*

Infection control kit, XO FLOW

- 2 * Instrument holder cover (MN-621)
- 6 * Tip for syringe (SD-214)
- 6 * Handle (AP-732)
- 1 * Cover for high-volume suction nozzle (AP-714)
- 1 * Cover for saliva suction nozzle (AP-715)
- 6 * Suction filter with holder (AP-963)
- 1 * Gold trap for cuspidor (MG-150)
- 1 * Cover for gold trap (MG-894)

Dashboard protection foil (to be delivered [from](#)).

- Protects the Dashboard from being scratched
- Self-adhesive foil covering the whole dashboard
- Cut-out for the instrument holder